

Halogen Moisture Analyzer

RLIFT-MCA110-10 RLIFT-MCA110-5 RLIFT-MCA110-2
RLIFT-MCA110-1 RLIFT-MCA110-1A



Description

- 

Halogen lamp/HBM sensor
- 

The percentage of moisture content/dry residual
- 

Aluminum shell/Stainless steel chamber
- 

Store historical set
- 

LCD(White back light with black font)
- 

Option:Aluminum box/Printer/Grinder
- 

Temperature/Time set

Product details



Specifications

Model	RLIFT - MCA110-10	RLIFT - MCA110-5	RLIFT - MCA110-2	RLIFT - MCA110-1	RLIFT - MCA110-1A
Max. capacity	110g				
Readability	0.01g	0.005g	0.002g	0.001g	
Operating Temp	5℃-35℃				
Display	LCD(White back light with black font)				LCD HD Back Display
Calibration	External Calibration				
Interface	RS232/RJ45/USB...				
Pan Size	Ø90mm				
Packing Size	490x350x360mm				
Gross Weight	7.5kg				
Heating Source	Halogen Lamp-1				Halogen Lamp-2
Temp.Sensor	PT100				
Temp.Range	40℃-199℃				
Temp.Step	1℃				
Moisture Range	0.00%-100.00%				
Moisture Readability	0.20%	0.10%	0.04%	0.01%	
Dry Range	100.00%-0.00%				
Dry Readability	0.20%	0.10%	0.04%	0.01%	
Temp Setting	40℃-199℃ By 1℃ Step				
Time Setting	1-99min By 10s				
Number Of Storage	15				
Heater	220V ±15% 50HZ/110V±15% 60HZ/400W				
Balance	Input:220V 50HZ/110V 60HZ Output:9V				