

Multifunctional Ultrasonic Cleaner

RLIFT-USC0140M RLIFT-USC0240M RLIFT-USC0340M
RLIFT-USC0540M RLIFT-USC0440M RLIFT-USC0640M
RLIFT-USC0740M RLIFT-USC1040M RLIFT-USC1540M
RLIFT-USC2040M RLIFT-USC2240M RLIFT-USC3040M
RLIFT-USC4540M



Description

- ◆ With sleeping mode to keep electrical energy saving
- ◆ Memory function will automatically save the most recently used setup parameters for next use.
- ◆ Ultra-low noise, adjustable time, adjustable power, frequency sweep function, degassing function, memory function, sleeping mode.
- ◆ Time setting with 1-99 minutes and built-in heating to 60℃, adjustable power range of 10%-100%.

Specifications

Model	RLIFT - USC0340M	RLIFT - USC0540M	RLIFT - USC0440M	RLIFT - USC0640M	RLIFT - USC0740M	RLIFT - USC1040M	RLIFT - USC1540M	RLIFT - USC2040M	RLIFT - USC2240M	RLIFT - USC3040M	RLIFT - USC4540M
Ultrasonic Power(W)-Max	120	120	180	180	240	240	360	360	480	600	720
Ultrasonic Power setting-%	10%-100%(step by 10%)										
Heating Power(W)-Max	100	100	200	200	500	500	500	500	500	500	1000
Capacity (L)-Max	3.2	4.8	4.5	6.5	7	10	15	20	22	30	45
Tank Size(L*W*H)-MM	240*135*100	240*135*150	300*150*100	300*150*150	300*240*100	300*240*150	330*300*150	330*300*200	500*300*150	500*300*200	500*300*300
External Size(L*W*H)-MM	265*165*235	265*165*285	330*180*260	330*180*310	330*270*260	330*270*310	360*330*330	360*330*380	530*330*340	530*330*390	530*330*490
Transducer Qty(PCS)	2	2	3	3	4	4	6	6	8	10	12
Drain Valve	Yes										
Carrying Handle	Yes										
Ultrasonic Freq(KHz)	40										
Temperature Setting(℃)-Max	Default 60℃, optional 80℃ (free of charge, please note when ordering)										
Timer Setting(MINS)	1-99										
Cleaning Basket	Yes										

Features



With large LCD display to show the precise time, temperature, and power level.



Humanized carrying handle and drain valve for whole model.



Built-in “sound insulation paste” keeps this series with low noise. With sweep function & degas function.

Specifications

Model	RLIFT - USC0140M	RLIFT - USC0240M
Ultrasonic Power(W)-Max	70	70
Ultrasonic Power setting-%	10%-100%(step by 10%)	
Heating Power(W)-Max	100	100
Capacity (L)-Max	1.3	2
Tank Size(L*W*H)-MM	150*135*65	150-135*100
External Size(L*W*H)-MM	180*165*205	180*165*205
Transducer Qty(PCS)	1	1
Drain Valve	No	
Carrying Handle	No	
Ultrasonic Freq(KHz)	40	
Temperature Setting(℃)-Max	60	
Timer Setting(MINS)	1-99	
Cleaning Basket	Yes	

Selling points

- ◆ Adjustable power
- ◆ Low noise
- ◆ Sweep function
- ◆ Degassing function
- ◆ Heating function
- ◆ Adjustable time
- ◆ Memory function
- ◆ Sleeping mode