

Water Bath

RLIFT-WB10-A RLIFT-WB30-A RLIFT-WB10-B RLIFT-WB30-B



- The chamber, frame and top cover are made of stainless steel.
- Between the inner and outer space are filled with high-quality insulation materials,anti-corrosive and good appearance.

Features



- Microprocessor control with timing and temperature correction function.
- LCD display
- The shell is made of high quality steel plate ,surface coated.

Specifications

Model	RLIFT-WB10-A	RLIFT-WB30-A	RLIFT-WB10-B(3 Usage)	RLIFT-WB30-B(3 Usage)
Temperature Range (℃)	RT+5~95	RT+5~95	RT+5~100	RT+5~100
Temperature Fluctuation (℃)	±0.5	±0.5	±0.5	±0.5
Temperature Uniformity(℃)	±0.5	±0.5	±0.5	±0.5
Power Supply	220V /50Hz	220V /50Hz	220V /50Hz	220V /50Hz
Power Rating (W)	1500±10%	2000±10%	1500±10%	2000±10%
Chamber Size (W×D×H)cm	30×24×20	50×30×20	30×24×20	50×30×20
Volume(L)	14	30	14	30
Package Size (W×D×H)cm	55×38×35	87×47×47	55×38×35	87×47×47
Net/Gross Weight(kg)	7.5/14	11.5/21	8.5/15	12.5/22