




# HYPODERMIC SYRINGE WITHOUT NEEDLE









## SPECIFICATIONS:

- ✓ Manufactured using high clarity barrel made from polypropylene.
- ✓ Proven for safety & free from microbes.
- ✓ Its latex-free piston is lubricated with medical grade silicone for a smooth draw.
- ✓ Internal ring in the barrel acts as a lock to ensure the plunger does not displace during draw.
- ✓ Exact long-lasting graduation mark printed on the barrel.
- ✓ Sterilized by ETO & its sterility is guaranteed
- ✓ Manufactured according to Standard ISO 7886-1 & ISO 11135



REGISTERED UNDER ACT 737  
MDA Registration No. : GA8289892816  
MDAMD 0061/2022

## DATASHEET:

Size	Type / Product Code			  Pyrogen Free  Latex Free  Single Use <div>STERILE EO</div> ETO Sterilized  Non-Toxic
	 Slip Center	 Luer Lock	 Slip Eccentric	
1ml	HSSC01	-	-	
3ml	HSSC03	HSLL03	-	
5ml	HSSC05	HSLL05	-	
10ml	HSSC10	HSLL10	HSSE10	
20ml	HSSC20	HSLL20	HSSE20	
50ml	HSSC50	HSLL50	HSSE50	

Distributed By:

**INTANIAGA SDN BHD** (658891-X)

Lot 2, Jalan CJ 1/4, Kawasan Perusahaan Cheras Jaya, 43200 Cheras Selangor, Malaysia.

Scan for more  
information

