



Centrifuge Tubes

ROYLAB Centrifuge Tubes are specially curated to fulfil lab users' needs for sample storage and centrifugation.

Product specifications:

- ▶ Made of high-quality polypropylene (PP)
- ▶ Leak-proof cap with large ridges for optimal grip
- ▶ ¾ turn closure for one-handed operation
- ▶ Autoclavable (121°C, 20 minutes); withstand temperature from -80°C to 120°C; tested to endure up to 16,000 RCF
- ▶ Clear, easy-to-read graduation marks & a large writing area
- ▶ Pack size: 25 pieces/ bag
- ▶ Volume: 15 or 50 ml
- ▶ DNase-free, RNase-free and non-pyrogenic

Ordering Information

Cat No.	Product Description	Pack Size
GLC-CT15	RoyLab 15ml Conical Centrifuge Tube	1 bag = 25 pieces; 1 carton = 20 bags
GLC-CT50	RoyLab 50ml Conical Centrifuge Tube	1 bag = 25 pieces; 1 carton = 20 bags



RoyLab Global Pte Ltd



enquiry@roylabglobal.com.sg



+65 8487 7719