

Product features:

Universal fit

- Compatible with virtually any single, multi-channel or electronic pipette brands on the market

Clean claims

- Pipette tips are certified free of DNase, RNase & pyrogen contamination

No presence of demolding aids, additives like DiHEMA & Oleamid are detected that might potentially interfere with biological tests, rendering false results being reported



Maximum sample recovery

- Ultra-low retention tips with low surface tension reduce sample hang-up, improve pipetting accuracy & minimise sample loss significantly when dealing with viscous media



Ultra Low Retention Tip



Brand X Standard Tip



Superior quality

- Stringent quality control procedures to ensure lot-to-lot consistency
- ISO 9001:2008 & CE certified (Each batch of tips is furnished with Certificate of Analysis)
- High-density filter prevents aerosol & liquid contamination
- Long reach design to prevent cross-contamination
- Precision injection moulding care for an airtight fit, which is critical for precision pipetting



Varying packaging format

- Various volume range covering from 10–1250µl
- Different packaging styles: Sterile, non-sterile, filter, non-filter, extended length, wide bore or bevel-shaped
- Packing available in racks, bulk, refills (with colour-coded tip platforms for easy identification)
- Graduated tips are meant for quick volume check

*Non-filter tips are autoclavable (121°C, 15 mins, 15 psi) for sterilisation purposes



RoyLab pipette tips are precision engineered to provide accurate and consistent results. With a secure, leak-proof fit and a range of sizes to choose from, our pipette tips are the essential choice for professionals who demand the best.



RoyLab Global Pte Ltd
enquiry@roylabglobal.com.sg
+65 8487 7719



Pipette
Tips

ORDERING INFORMATION

Catalogue No.	Series	Product Specifications						Unit	Packaging	
		Colour	Volume	Retention	Filter	Sterile	Style		Pack Size	Description
P-PIP-S-NF-RAC-S-10-SS	Aqua	Natural	10µl	Standard	No	Yes	Short shank	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Sterile, Standard Retention, Non-filtered, Short Shank, Racked
P-PIP-S-NF-RAC-S-10-LS		Natural	10µl	Standard	No	Yes	Long shank	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Sterile, Non-filtered, Long Shank, Racked
P-PIP-S-NF-RAC-S-200		Natural	200µl	Standard	No	Yes	Long shank	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Sterile, Non-filtered, Racked
P-PIP-S-NF-RAC-S-300		Natural	300µl	Standard	No	Yes	NA	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Sterile, Non-filtered, Racked
P-PIP-S-NF-RAC-S1250-S		Natural	1250µl	Standard	No	Yes	NA	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Sterile, Non-filtered, Racked
P-PIP-S-NF-RL-S-10-SS		Natural	10µl	Standard	No	Yes	Short shank	Tray (Reloads)	"100(1 tray = 96 tips; 1 box = 10 trays; 1 carton = 10 boxes)"	Sterile, Non-filtered, Short Shank, Reloads
P-PIP-S-NF-RL-S-10-LS		Natural	10µl	Standard	No	Yes	Long shank	Tray (Reloads)	"100(1 tray = 96 tips; 1 box = 10 trays; 1 carton = 10 boxes)"	Sterile, Non-filtered, Long Shank, Reloads
P-PIP-S-NF-RL-S-200		Natural	200µl	Standard	No	Yes	NA	Tray (Reloads)	"100(1 tray = 96 tips; 1 box = 10 trays; 1 carton = 10 boxes)"	Sterile, Non-filtered, Reloads
P-PIP-S-NF-RL-S-300		Natural	300µl	Standard	No	Yes	NA	Tray (Reloads)	"100(1 tray = 96 tips; 1 box = 10 trays; 1 carton = 10 boxes)"	Sterile, Non-filtered, Reloads
P-PIP-S-NF-RL-S-1250-S		Natural	1250µl	Standard	No	Yes	NA	Tray (Reloads)	"80 (1 tray = 96 tips; 1 box = 8 trays; 1 carton = 10 boxes)"	Sterile, Non-filtered, Reloads
P-PIP-NS-NF-B-S-10-SS		Natural	10µl	Standard	No	No	Short shank	Bag	"10 (1 bag = 1000 tips; 1 carton = 10 bags)"	Non-Sterile, Non-filtered, Short Shank, Bagged
P-PIP-NS-NF-B-S-10-LS		Natural	10µl	Standard	No	No	Long shank	Bag	"10 (1 bag = 1000 tips; 1 carton = 10 bags)"	Non-Sterile, Non-filtered, Long Shank, Bagged
P-PIP-NS-NF-B-S-200		Natural	200µl	Standard	No	No	NA	Bag	"10 (1 bag = 1000 tips; 1 carton = 10 bags)"	Non-Sterile, Non-filtered, Bagged
P-PIP-NS-NF-B-S-300		Natural	300µl	Standard	No	No	NA	Bag	"10 (1 bag = 1000 tips; 1 carton = 10 bags)"	Non-Sterile, Non-filtered, Bagged
P-PIP-NS-NF-B-S-1250-S		Natural	1250µl	Standard	No	No	NA	Bag	"10 (1 bag = 1000 tips; 1 carton = 10 bags)"	Non-Sterile, Non-filtered, Bagged
P-PIP-S-F-S-10-SS		Natural	10µl	Standard	Yes	Yes	Short shank	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Sterile, Filtered, Short Shank, Racked
P-PIP-S-F-S-10-LS	Glacier	Natural	10µl	Standard	Yes	Yes	Long shank	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Sterile, Filtered, Long Shank, Racked
P-PIP-S-F-S-20		Natural	20µl	Standard	Yes	Yes	NA	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Sterile, Filtered, Racked
P-PIP-S-F-S-100		Natural	100µl	Standard	Yes	Yes	NA	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Sterile, Filtered, Racked
P-PIP-S-F-S-200		Natural	200µl	Standard	Yes	Yes	NA	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Sterile, Filtered, Racked
P-PIP-S-F-S-1250-S		Natural	1250µl	Standard	Yes	Yes	NA	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Sterile, Filtered, Racked
P-PIP-S-NF-RAC-P-10-SS	Icicle	Natural	10µl	Ultra-low	No	Yes	Short shank	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Premium, Ultra Low Retention, Sterile, Non-filtered, Short Shank, Racked
P-PIP-S-NF-RAC-P-10-LS		Natural	10µl	Ultra-low	No	Yes	Long shank	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Premium, Ultra Low Retention, Sterile, Non-filtered, Long Shank, Racked
P-PIP-S-NF-RAC-P-200		Natural	200µl	Ultra-low	No	Yes	NA	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Premium, Ultra Low Retention, Sterile, Non-filtered, Racked
P-PIP-S-NF-RAC-P-300		Natural	300µl	Ultra-low	No	Yes	NA	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Premium, Ultra Low Retention, Sterile, Non-filtered, Racked
P-PIP-S-NF-RAC-P-1250-S		Natural	1250µl	Ultra-low	No	Yes	NA	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Premium, Ultra Low Retention, Sterile, Non-filtered, Racked
P-PIP-S-NF-RL-P-10-SS		Natural	10µl	Ultra-low	No	Yes	Short shank	Tray (Reloads)	"100 (1 tray = 96 tips; 1 box = 10 trays; 1 carton = 10 boxes)"	Premium, Ultra Low Retention, Sterile, Non-filtered, Short Shank, Reloads
P-PIP-S-NF-RL-P-10-LS		Natural	10µl	Ultra-low	No	Yes	Long shank	Tray (Reloads)	"100 (1 tray = 96 tips; 1 box = 10 trays; 1 carton = 10 boxes)"	Premium, Ultra Low Retention, Sterile, Non-filtered, Long Shank, Reloads
P-PIP-S-NF-RL-P-200		Natural	200µl	Ultra-low	No	Yes	NA	Tray (Reloads)	"100 (1 tray = 96 tips; 1 box = 10 trays; 1 carton = 10 boxes)"	Premium, Ultra Low Retention, Sterile, Non-filtered, Reloads
P-PIP-S-NF-RL-P-300		Natural	300µl	Ultra-low	No	Yes	NA	Tray (Reloads)	"100 (1 tray = 96 tips; 1 box = 10 trays; 1 carton = 10 boxes)"	Premium, Ultra Low Retention, Sterile, Non-filtered, Reloads
P-PIP-S-NF-RL-P-1250-S		Natural	1250µl	Ultra-low	No	Yes	NA	Tray (Reloads)	"80 (1 tray = 96 tips; 1 box = 8 trays; 1 carton = 10 boxes)"	Premium, Ultra Low Retention, Sterile, Non-filtered, Reloads
P-PIP-NS-NF-B-P-10-SS		Natural	10µl	Ultra-low	No	No	Short shank	Bag	"10 (1 bag = 1000 tips; 1 carton = 10 bags)"	Premium, Ultra Low Retention, Non-Sterile, Non-filtered, Short Shank, Bagged
P-PIP-NS-NF-B-P-10-LS		Natural	10µl	Ultra-low	No	No	Long shank	Bag	"10 (1 bag = 1000 tips; 1 carton = 10 bags)"	Premium, Ultra Low Retention, Non-Sterile, Non-filtered, Long Shank, Bagged
P-PIP-NS-NF-B-P-200		Natural	200µl	Ultra-low	No	No	NA	Bag	"10 (1 bag = 1000 tips; 1 carton = 10 bags)"	Premium, Ultra Low Retention, Non-Sterile, Non-filtered, Bagged
P-PIP-NS-NF-B-P-300		Natural	300µl	Ultra-low	No	No	NA	Bag	"10 (1 bag = 1000 tips; 1 carton = 10 bags)"	Premium, Ultra Low Retention, Non-Sterile, Non-filtered, Bagged
P-PIP-NS-NF-B-P 1250-S		Natural	1250µl	Ultra-low	No	No	NA	Bag	"10 (1 bag = 1000 tips; 1 carton = 10 bags)"	Premium, Ultra Low Retention, Non-Sterile, Non-filtered, Bagged
P-PIP-S-F-P-10-SS	Poseidon	Natural	10µl	Ultra-low	Yes	Yes	Short shank	Rack	"50 (1 rack=96 tips; 1 box=10 racks; 1 carton=5 boxes)"	Premium, Ultra Low Retention, Sterile, Filtered, Short Shank, Racked
P-PIP-S-F-P-10-LS		Natural	10µl	Ultra-low	Yes	Yes	Long shank	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Premium, Ultra Low Retention, Sterile, Filtered, Long Shank, Racked
P-PIP-S-F-P-20		Natural	20µl	Ultra-low	Yes	Yes	NA	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Premium, Ultra Low Retention, Sterile, Filtered, Racked
P-PIP-S-F-P-100		Natural	100µl	Ultra-low	Yes	Yes	NA	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Premium, Ultra Low Retention, Sterile, Filtered, Racked
P-PIP-S-F-P-200		Natural	200µl	Ultra-low	Yes	Yes	NA	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Premium, Ultra Low Retention, Sterile, Filtered, Racked
P-PIP-S-F-P-300		Natural	300µl	Ultra-low	Yes	Yes	NA	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Premium, Ultra Low Retention, Sterile, Filtered, Racked
P-PIP-S-F-P-1250-S		Natural	1250µl	Ultra-low	Yes	Yes	NA	Rack	"50 (1 rack = 96 tips; 1 box = 10 racks; 1 carton = 5 boxes)"	Premium, Ultra Low Retention, Sterile, Filtered, Racked
P-PIP-R10	NA	Natural	Empty rack for 10µl pipette tips					Rack	"50 (1 box = 10 racks; 1 carton = 5 boxes)"	Empty Rack for 10ul
P-PIP-R200		Natural	Empty rack for 200µl pipette tips					Rack	"50 (1 box = 10 racks; 1 carton = 5 boxes)"	Empty Rack for 200ul
P-PIP-R300		Natural	Empty rack for 300µl pipette tips					Rack	"50 (1 box = 10 racks; 1 carton = 5 boxes)"	Empty Rack for 300ul
P-PIP-R1250		Natural	Empty rack for 1250µl pipette tips					Rack	"50 (1 box = 10 racks; 1 carton = 5 boxes)"	Empty Rack for 1250ul

[illegible]